

## Post and Core Restoration System

This complete post and core restoration system is designed to provide optimal restorations. The system includes LuxaCore Z-Dual post cement, for exceptional compressive strength, cuttability and polishability similar to dentin as well as optimal flow; LuxaBond–Total Etch dual-curing adhesive, suited for indirect techniques; and LuxaPost composite post, which reacts to stress like dentin.



**Zenith Dental**  
800-662-6383  
[www.zenithdental.com](http://www.zenithdental.com)  
**Select 10.**

## Veneer Delivery System

Featuring a proprietary tray delivery system for individually made hybrid or pure porcelain veneers, the GlamSmile veneer system is designed to make time-consuming veneer placement more efficient. The pain-free system does not require removal of healthy tooth structure. Other features include minimal to no prep, precise design, and minimal cost when compared with traditional veneers. The digital design system can create a full arch (up to



10 teeth) per 1-hour application.  
**Remedent Inc.**  
888-913-9899  
[www.glamsmile.com](http://www.glamsmile.com)  
**Select 14.**

## Reception Area Educational DVD



CAESY Smile Channel DVD 2008 is the latest version of the educational reception area program designed to inform and entertain patients on a variety of topics related to general and esthetic dentistry. It now features two new programs: *Talk Out Loud* and *Health Watch*, which both

motivate patients to take a proactive role in maintaining their oral health. *Talk Out Loud* details important oral health topics for teens—such as smoking and nutrition—to initiate conversation between young patients and their dentists and hygienists. *Health Watch*, an interview-style news program, presents oral health as a component of total body health and encourages patients to request additional information and engage in conversations with their dental professionals. The programs feature 15 new segments that incorporate before-and-after images.

**Patterson Dental Supply Inc.**  
800-294-8504  
[www.caesy.com](http://www.caesy.com)  
**Select 11.**

## Composite System

The Filtek LS Low Shrink Posterior Restorative System is formulated to achieve less than 1% polymerization shrinkage and reportedly offers optimal handling, strong mechanical properties, and natural-looking restorations. Key to the composite system is silorane chemistry, which uses ring-opening polymerization to minimize shrinking. In addition, the silorane technology increases ambient light stability with up to 9 minutes of handling time under operatory light. The system also promises easy placement,



even in horizontal layers, and sculptability to facilitate the creation of occlusal anatomy before cure.

**3M ESPE**  
800-634-2249  
[www.3MESPE.com](http://www.3MESPE.com)  
**Select 12.**

## Handpiece Lubrication System


The ClearView automatic handpiece lubrication system, which features a compact design, is said to lubricate in 10 seconds with no overspray or odor. A built-in reservoir helps ensure a precise amount of lubricant is delivered at the correct air pressure. Oil and patient debris are captured to protect staff and eliminate counter cleanup.

The lubricator features one moving part, making it economical and reliable. In addition, it is air-powered, lubricates all air-driven handpieces, and comes with a 2-year warranty.

**Dental Air Solutions**  
877-212-9411  
[www.dentalairsolutions.com](http://www.dentalairsolutions.com)  
**Select 13.**

## Perio Probe

Used with the DIAGNOdent caries detection device, the DIAGNOdent perio probe is designed to detect calculus concretions in periodontal pockets up to 9 mm deep. The probe uses laser fluorescence of subgingival calculus to determine quantification of calculus left behind by standard root planning and scaling. Because calculus fluoresces differently than healthy tissue, the probe detects this difference, sends a visual signal to



its numerical display, and also offers an audible signal to indicate when calculus is identified. This increases patient awareness of their calculus levels, leading to increased compliance with recommended treatment. The device is available for use with DIAGNOdent Classic and pen versions.

**KaVo Dental**  
888-ASK-KAVO (888-275-5286)  
[www.kavousa.com](http://www.kavousa.com)  
**Select 15.**

## USB Sensor Controller

Designed to offer choice in sensor connectivity, the DEXusb USB sensor controller offers a fast, robust option for digital x-ray imaging. Building on the ultra-portability concept of its PC card, the controller allows clinicians to use a PC card or USB. The unit uses a mini-B USB connector, the de facto standard on portable devices such as digital cameras and cell phones. Designed for use with the DEXIS Perfect-Size sensor, the controller's lightweight, ergonomic design allows easy transport between operatories. The unit is self-pow-



ered and presents an economical way to bring digital x-ray to a practice. It can be used on a countertop or wall mounted.

**Dexis LLC**  
888-883-3947  
[www.dexis.com](http://www.dexis.com)  
**Select 16.**

## Fluoride-Releasing Bonding System



FL-Bond II fluoride-releasing bonding system is formulated with a combination primer/bonding agent that provides an exceptional bond to enamel and dentin alike with a secure marginal seal. Made with second-generation Giomer technology, FL-Bond II is regarded as durable

and reinforces the bonding interface between the restorative material and tooth structure. The HEMA-free primer prevents the deterioration of the bonding interface due to hydrolysis and ensures long-term durability of the bond. Easy to handle, the agent's viscosity ensures uniform application. The system's high level

of fluoride release and its anti-bacterial effect reportedly help protect against plaque formation.  
**Shofu**  
800-827-4638  
[www.shofu.com](http://www.shofu.com)  
**Select 17.**